

## MagicLink™ HRP Antibody Labeling Kits (100 µg, no desalt)

### Components


Components		Product size		Storage condition
		BP-50002	BP-50001	
		1 x 100 µg	3 x 100 µg	
A	MAGIC NHS (MW ~900)	1 vial	3 vials	-20C
B	LINK activated HRP	1 vial	3 vials	-20C
C	Reaction Buffer	10 ml	30 ml	RT

Note: The kit above is designed for IgG antibodies, but works well for any amine containing biomolecule. Please follow same technical tips if required.

### Overview

MagicLink™ HRP Antibody Labeling Kits are the 3rd generation HRP conjugation technology which can be used to conjugate horseradish peroxidase (HRP) to protein, antibody, amine modified oligo, etc. The labeling kits feature the most stable linkage between HRP antibody on the market. The instant and efficient labeling reactions yield 95 to 100% HRP conjugates.

These kits are specifically optimized to label antibodies at a scale for 100 ug. The kits' format is based on instant reaction between functional group MAGIC and LINK at room temperature. By following the protocol provided in the kits, the end users can label their antibody with MAGIC NHS to get MAGIC-antibody which instantly reacts with LINK activated HRP (provided in the kits) to achieve Ab-HRP conjugates.



## At a Glance

### Protocol summary

1. Add reaction buffer into antibody.
2. Transfer whole antibody solution to MAGIC NHS vial.
3. Incubate at room temperature for 60 minutes.
4. Mix MAGIC-antibody with LINK activated HRP.
5. Storage.

Note: Upon receipt, store the kit at 4°C. When stored properly, the kit should be stable for six months. Alternatively, components A and C can be stored at -20°C. Do not freeze reaction buffer (Component C). Warm all the components and centrifuge the vials briefly before opening, and immediately prepare the required solutions before starting your conjugation. The following SOP is an example for labeling goat anti-mouse IgG antibody.

## Preparation of Working Solution

For labeling 100 µg antibody (assuming the target antibody concentration is 2 mg/mL), final volume would be in 50 µl reaction buffer (component C), or 1X PBS pH 7.2 – 7.5. For lyophilized antibody, add 50 µl of reaction buffer to reconstitute. Antibody, in liquid form and not in 1X PBS, pH 7.2 – 7.5 may require buffer exchange into reaction buffer, see note.

### Note

- If you have a different concentration, adjust the antibody volume accordingly to make ~100 µg antibody available for your labeling reaction.
- The antibody should be dissolved in 1X phosphate buffered saline (PBS), pH 7.5; if the antibody is dissolved in glycine buffer, it must be dialyzed against 1X PBS, pH 7.2-7.5. Alternatively, use Amicon Ultra-0.5, Ultracel-10 Membrane, 10 kDa (Cat # UFC501008 from Millipore) to remove free amines or ammonium salts (such as ammonium sulfate and ammonium acetate) that are widely used for antibody precipitation, and buffer exchange into reaction buffer.
- Impure antibodies or antibodies stabilized with bovine serum albumin (BSA) or gelatin will not be labeled well.
- The conjugation efficiency is significantly reduced if the antibody concentration is less than 1 mg/ml. For optimal labeling efficiency the final antibody concentration range of 1-5 mg/ml is recommended.
- The presence of sodium azide will inhibit HRP activity.

## Labeling Protocol

React Antibody with MAGIC NHS reaction:

Add the antibody solution directly into the vial of MAGIC NHS (Component A), and mix them well by repeatedly pipetting for a few times or vortex the vial for a few seconds. Keep the antibody labeling reaction mixture at room temperature for 60 minutes. The antibody-labeling reaction mixture can be rotated or shaken for longer time if needed.

## HRP-Antibody Conjugation

1. Make LINK-activated HRP solution by adding 50  $\mu$ L ddH<sub>2</sub>O into the vial of HRP-LINK (Component B), mix well by repeatedly pipetting for a few times or vortex the vial for a few seconds.
2. Mix LINK-HRP solution with MAGIC-antibody solution at different ratio, mix well and rotating the mixture for 1 hour at room temperature.

### Note:

- It is recommended the mix link activated HRP and antibody at 3:1.
- Use all 50  $\mu$ l of LINK-HRP to label 100  $\mu$ g MAGIC-antibody (generally 150 kDa) at 5x HRP to antibody ratio mole ratio, 40  $\mu$ L for 4x, 20  $\mu$ L for 2x, etc.
- Magic Activated protein/antibody should be used right away.
- For a different protein, user need optimize the Protein/HRP mix ratio for fit the application accordingly.

The HRP-antibody conjugate is now ready to use. For immediate use, the HRP-antibody conjugate need be diluted with the buffer of your choice. For longer term storage, HRP-antibody conjugate solution need be concentrated or freeze dried.

## Storage of HRP-Antibody Conjugate

The antibody conjugate should be stored at > 0.5 mg/mL in the presence of a carrier antibody (e.g., 0.1% bovine serum albumin). For longer storage, the HRP-antibody conjugates could be lyophilized and stored at  $\leq -20$  °C.

## Troubleshooting

Problem	Possible cause	solution
Low or no conjugation with MAGIC NHS	Buffer containing primary amine	Buffer exchange the antibody into a non-amine-containing buffer such as the PBS provided by desalting columns or dialysis
	MAGIC NHS was hydrolyzed	Use reagent immediately upon reconstitution
	Carrier protein was present in the antibody solution	Remove carrier protein before conjugation by using Protein A, G or A/G resin or an antibody clean-up kit. This will reduce competition for labeling