

MagicLink™ Biotin Antibody Labeling Kit: Sulfo-NHS-LC-Biotin

Components


Component		Product size		Storage
		BP-50058	BP-50059	
		8 rxn	8 rxn	
A	Sulfo-NHS-LC-Biotin	8 x 2 mg	8 x 2 mg	-20C
B	Reaction buffer	30 ml	N/A	4-8°C
C	Protein concentrator	8	N/A	RT

Note: BP-50059, user may use 1x PBS pH 7.2 – 7.4 as reaction buffer is not supplied.

Overview

BroadPharm Sulfo-NHS-LC-biotin Antibody Biotinylation Kit provides optimized reagents for labeling antibodies and protein concentrator for purifying the labeled molecule. Each reaction is sufficient for labeling 50-200 µg of antibody in 100 µl reaction volumes.

The Sulfo-NHS-LC-Biotin is an intermediate-length, water-soluble biotinylation reagent for labeling antibodies, proteins and other molecules that have primary amines. This kit is specifically optimized to label antibodies at a scale up to 1 mg. The kit format is a convenient single-use microtubes, eliminating difficulties associated with weighing small quantities of reagent. Biotin is a small, naturally occurring vitamin that binds with high affinity to avidin and avidin-like proteins. Biotinylated antibodies typically retain biological activity because the biotin group is relatively small. An antibody conjugated with several biotin molecules can interact rapidly and tightly with streptavidin. N-Hydroxysuccinimide (NHS) esters are the most popular biotinylation reagents. In pH 7-9 buffers, NHS esters react efficiently with primary amino groups (-NH₂) by nucleophilic attack, forming an amide bond and releasing the NHS. Proteins typically have many sites for labeling, including the primary amine in the side chain of lysine (K) residues and the N-terminus of each polypeptide.



Sulfo-NHS-LC-Biotin
Mol. Wt.: 556.59
Spacer Arm: 22.4A

At a Glance

Protocol summary

1. Add 200 μ l DI H₂O to a Sulfo-NHS-LC-Biotin vial.
2. Prepare the antibody at 2 mg/ml in reaction buffer or pH 7.4 PBS buffer.
3. Add the calculated amount of biotin working solution to antibody stock solution.
4. Incubate at room temperature for 30 - 60 minutes.
5. Remove excess biotin reagent with protein concentrator.

Note: Upon receipt, store the reaction buffer at 4°C, and component Sulfo-NHS-LC-Biotin at -20°C. When stored properly, the kit should be stable for six months.

Do not freeze component B reaction buffer. Warm all the components to room temperature and centrifuge the vials briefly before opening, and immediately prepare the required solutions before starting your conjugation. The following SOP is an example for labeling goat anti-mouse IgG antibody.

Preparation of Working Solution

1. Antibody working solution:

For labeling 1 mg antibody, the preferred concentration antibody concentration is 2 mg/ml.

Note:

- If you have a different concentration, adjust the antibody concentration accordingly
- The antibody should be dissolved the reaction buffer. If the antibody is dissolved in glycine buffer, it must be dialyzed against 1X PBS, pH 7.2-7.4, or use Amicon Ultra-0.5, Ultracel-10 Membrane, 10K MWCO (Cat # UFC501008 from Millipore) to remove free amines or ammonium salts (such as ammonium sulfate and ammonium acetate) that are widely used for antibody precipitation.
- Impure antibodies or antibodies stabilized with bovine serum albumin (BSA) or gelatin will not be labeled well.
- The conjugation efficiency is significantly reduced if the antibody concentration is less than 1 mg/ml. For optimal labeling efficiency the final antibody concentration range of 1-5 mg/ml is recommended.

2. Sulfo-NHS-LC-biotin working solution

Add 200 μ L Di-H₂O into the one vial of Sulfo-NHS-LC-Biotin, the concentration of Sulfo-NHS-LC-Biotin is 2mg/200ul.

Note: Keep the Sulfo-NHS-LC-Biotin working solution in ice/H₂O, and use it as soon as possible. Unused stock solution should be discarded after the conjugation.

3. Calculation

The Degree of Labeling (DOL) depends on the size and distribution of amino groups on the protein and the amount of biotin reagent used. Compared to reactions involving concentrated protein solutions, labeling reactions with dilute protein solutions require a greater fold molar excess of biotin reagent to achieve the same incorporation level. Experiments that used a 20-fold molar excess of biotin reagent to label 1-5 mg/ml antibody (IgG) resulted in 4-6 biotin groups per antibody molecule. Experiments that used a 50-fold molar excess of biotin reagent to label 50-200 μ g of antibody (in 200-700 μ l) resulted in 1-3 biotin groups per antibody molecule. Adjust the molar ratio of Sulfo-NHS-LC-Biotin to protein to obtain the desired level of incorporation.

Formula:

$$\begin{aligned}
 V_{\text{biotin_NHS}}(\mu\text{l}) &= \frac{\text{Amount_biotin_NHS}(\text{mmol}) \times \text{MW_biotin_NHS}(\text{mg}/\text{mmol})}{\text{Conc_biotin_NHS}(\text{mg}/\mu\text{l})} \\
 &= 20 \times \frac{\text{Amount_protein}(\text{mmol}) \times \text{MW_biotin_NHS}(\text{mg}/\text{mmol})}{\text{Conc_biotin_NHS}(\text{mg}/\mu\text{l})} \\
 &= 20 \times \frac{V_{\text{protein}}(\text{ml}) \times \text{Conc_protein}(\text{mg}/\text{ml})}{\text{MW_protein}(\text{mg}/\text{mmol})} \times \frac{\text{MW_biotin_NHS}(\text{mg}/\text{mmol})}{\text{Conc_biotin_NHS}(\text{mg}/\mu\text{l})} \\
 &= 20 \times \frac{V_{\text{protein}}(\text{ml}) \times \text{Conc_protein}(\text{mg}/\text{ml})}{\text{MW_protein}(\text{mg}/\text{mmol})} \times \frac{556(\text{mg}/\text{mmol})}{2/200(\text{mg}/\mu\text{l})}
 \end{aligned}$$

Example:

For 1 ml of a 2 mg/ml IgG (150,000 MW) solution, 14.8 μl of Sulfo-NHS-LC-Biotin will be added.

$$\begin{aligned}
 V_{\text{biotin_NHS}}(\mu\text{l}) &= 20 * 1 * 2 / 150000 * 556 / (2 / 200) \\
 &= 14.8 \mu\text{l}
 \end{aligned}$$

Conjugation Experimental Protocol

1. Add the amount of Sulfo-NHS-LC-Biotin needed (as calculated in Section Calculation) to the antibody solution.
2. Mix by gently pipetting up and down.
3. Incubate the reaction at room temperature for 30 minutes.

Excess Biotin Removal

1. Hydrate concentrator membrane 'filter device' with 400 to 500 μl of reaction buffer or DI water, and microcentrifuge 14,000 x g, for 3 minutes. Discard, liquid from filter device and collection tube.
2. Spin down by adding biotin labeled antibody to the concentrator/filter device up to 500 μl . Microcentrifuge at 14,000 x g, 8 minutes, or to minimum volume \sim 50 μl left in the filter device. Discard waste from the collection tube.
3. Desalt by adding reaction buffer to the filter device up to 500 μl . Microcentrifuge at 14,000 x g, 8 minutes, or to minimum volume \sim 50 μl left in the filter device. Discard waste from the collection tube.
4. Repeat step 3, twice.
5. Collect labeled antibody from filter device into a microcentrifuge tube.
6. Optional for maximum recovery, add a small volume of the reaction buffer (volume determined by the user), to the filter device to rinse out residual antibody, microcentrifuge pulse spin, collect antibody/reaction buffer from filter device, and add to the microcentrifuge tube from step 5, mix.
7. Store biotinylated antibody at 4°C for < 1 month. For longer periods, store at -20°C or -80°C (e.g., 0.1% bovine serum albumin) and 0.02-0.05% sodium azide.

Biotin Quantitation

Use BroadPharm's Biotin quantitation kit to calculate the degree of labeling (DOL), cat# BP-50060.

Storage of Antibody-Biotin

The antibody conjugate should be stored in the presence of a carrier protein.

Troubleshooting

Problem	Possible cause	Solution
Low or no biotinylation	Buffer containing primary amine	Buffer exchange the antibody into a non-amine-containing buffer such as the PBS provided by desalting columns or dialysis
	Sulfo-NHS-LC biotin was hydrolyzed	Use reagent immediately upon reconstitution
	Carrier protein was present in the antibody solution	Remove carrier protein before biotinylation by using Protein A, G or A/G resin or an antibody clean-up kit. This will reduce competition for labeling